

NEW

MADE BY TETRIX PatentDispenser

One dispenser for all pump bottles


+358 10 219 35 60

© TETRIX 2020

**URNS EVERYDAY PUMP
BOTTLES INTO PRO
LEVEL DISPENSERS**



Non-brand-loyal dispenser

PatentDispenser works with ANY supplier's soap, cream disinfection liquid, etc.

All sizes of bottles can easily be installed into the mechanism, solving the logistical problem of dispenser liquids.

Back up system.

Institutions often risk getting locked to purchase a certain manufacturer's products. If deliveries fail, the whole institutions stands without e.g. disinfection liquid. A back-up of PatentDispensers can overcome this dilemma in times of crisis.

Everyday bottles with a pump.

You can insert all kinds of bottles that have a pump, into the dispenser magazine. Easily adjust the level of the brackets to fit any pump bottle.

Don't throw away the bottles!

In case of a pandemic, even the bottles we consider to be "everyday stuff", might run out. Save them for use in your PatentDispenser!



SAVE MONEY - RECYCLE BOTTLES

Refill your old bottles and buy money saving container packs. We guarantee that you will save money, and be kinder to the environment!



Easy as a stapler

Insert your bottle into the magazine and adjust the dispenser in seconds. The smartest solutions are sometimes the most obvious ones. The dispenser is ideal for bottles 60-75mm in diameter, also e.g. oval shaped bottles. It is easy!

DesiStand & PatentDispenser.

Sanitizers located close to entries actually get used. Combine your PatentDispenser with DesiStand. Avoid the main obstacles and achieve better sanitation:

A) **PatentDispenser solves obstacles on a macro level.** When hitting e.g. a pandemic, there is a risk that producers run out of one-brand-dispenser systems or containers for them. The solution is PatentDispenser - as you can use any bottle with a pump and refill it.

B) **PatentDispenser solves obstacles on a micro level.** The risk is that people do not use e.g. disinfectants. The solution is a combo of PatentDispenser and DesiStand. This makes the sanitizer mobile and available in places where people will actually use it.

Don't invite **GERMS**

Bottles standing on a desk or a sink might be a good idea. Yet in crowded environments the bottle might just be an attraction for contaminated hands to touch the surfaces - making matters worse!

Use our dispenser to avoid unnecessary contamination!



Higher level of hygiene with stand- and wall mounting

Inserted into the PatentDispenser, your everyday bottle turns into a hospital-grade solution.

- Can be re-adjusted to any size in seconds
- Only highest-quality parts and finishing
- Cheaper refilling with bulk containers
- You are not locked to a choice of one distributor
- Robot precision of 0,02mm
- Electroplated surfaces
- Long elbow trigger for better hygiene
- Only 2 screws for rapid wall mounting
- Can be sterilized in a washing machine
- Works as back-up or standard dispenser
- Extras such as hangers and locking
- We can pre-adjust the dispenser before delivery

The PatentDispenser way.

If not now, what about the next time? The PatentDispenser works also fine between critical periods, making pump bottles hygienic. However, when the next peak of an outbreak comes (and it will) your PatentDispenser will save you many worries. Stay prepared! It's easy!



Order yours today!

The production is at full speed and we have the capacity to serve. If you think that our pump would be a solution for your hospital or other facility, feel free to contact us directly. We are also open to discuss with retailers that can make the PatentDispenser more accessible on their markets.

FLAT PACKED - SPACE SAVING

Thanks to the folding mechanism, the PatentDispenser is only about 3cm high when packed. A logistical ideal for customers further from us. Sold in packs of 10 pcs. The box is only about 30x30x10cm in size. Every parcel is Ethanol sprayed before closing, for added security against contamination.



The PatentDispensers are packed in boxes of 10 pcs. At present this is our minimum quantity order size - as we try to get the product out to as many users as possible. Please allow for hospitals and healthcare institutions to be served first. Production capacity at normal levels is 1000 units per day and can be raised according to need.

patentdispenser@tetrix.fi | +358 10 219 3560 | www.tetrix.fi